



# Hazardous Materials Management Sharp Containers Assembly & Packaging Guide

The following guide provides a brief overview regarding the assembly and disposal of 5-quart sharp containers provided by Hazardous Materials Management.





## Assembly and Use

Press the lid into the base making sure all six tabs are fully seated into the base (three in the front and three in the back) and have snapped fully into the base. The six tabs should be seen completely seated through the underside of the red base (pictured left). Please do not fill above the fill line (pictured left) as this may hinder the closure of the rotating lid and allow for sharp objects to protrude through.

## Disposal

- Submit the sharp container at <https://utsa.campusoptics.com/>, using the sharp container waste stream.
- Verify “rotating cylinder lid” tabs (2) are snapped into the lid.
- Place a piece of tape over the lid to further latch the sharp container (pictured below)
- Write the container ID# from Campus optics on the sharp container (pictured below).



## Notes

- Sharp containers are given out on a one-for-one basis
- No liquid-containing, ampules, test tubes, syringes, etc should be placed in the sharp containers. All liquid should be collected in another container and submitted as chemical waste.
- Failure to follow disposal procedures may result in rejection of waste submission.